

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat CDK1/2/3 Polyclonal Antibody

多克隆抗体

产品货号: PA02637

---

第三版

**描述:** Plays a key role in the control of the eukaryotic cell cycle. It is required in higher cells for entry into S-phase and mitosis. p34 is a component of the kinase complex that phosphorylates the repetitive C-terminus of RNA polymerase II.

**抗原:** Synthesized peptide derived from human Cdk1/2/3 around the non-phosphorylation site of Thr14.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** breast cancer tissues, cancer cell lines and glioblastoma tissue.