

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TRAIL Polyclonal Antibody

多克隆抗体

产品货号: PA02638

第三版

描述: Cytokine that binds to TNFRSF10A/TRAILR1, TNFRSF10B/TRAILR2, TNFRSF10C/TRAILR3, TNFRSF10D/TRAILR4 and possibly also to TNFRSF11B/OPG. Induces apoptosis. Its activity may be modulated by binding to the decoy receptors TNFRSF10C/TRAILR3, TNFRSF10D/TRAILR4 and TNFRSF11B/OPG that cannot induce apoptosis.

抗原: Synthesized peptide derived from the Internal region of human TRAIL.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widespread; most predominant in spleen, lung and prostate.