

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PRL Polyclonal Antibody

多克隆抗体

产品货号: PA02655

第三版

描述: Various hormones are secreted from the anterior pituitary during development and growth. Prolactin is a growth factor secreted by the anterior pituitary that is necessary for the proliferation and differentiation of the mammary glands.

抗原: Synthesized peptide derived from the Internal region of human Prolactin.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: