

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CUX1 Polyclonal Antibody

多克隆抗体

产品货号: PA02659

第三版

描述:The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified. CUX1 (Cut Like Homeobox 1) is a Protein Coding gene. Diseases associated with CUX1 include Fiedler's Myocarditis and Dirofilaria. Among its related pathways are Intra-Golgi traffic and FGFR1 mutant receptor activation. GO annotations related to this gene include sequence-specific DNA binding and RNA polymerase II regulatory region sequence-specific DNA binding. An important paralog of this gene is CUX2.

抗原:Synthesized peptide derived from the N-terminal region of human CUX1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: