
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CNTROB Polyclonal Antibody

多克隆抗体

产品货号: PA02667

第三版

描述: This gene encodes a centrosomal protein that interacts with BRCA2, and is required for centriole duplication and cytokinesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

抗原: Synthesized peptide derived from the Internal region of human Centrobilin

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed (at protein level). Highly expressed in testis. Also expressed in spleen, thymus, prostate, small intestine, colon and peripheral blood leukocytes.