

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RSAD2 Polyclonal Antibody

多克隆抗体

产品货号: PA02674

第三版

描述:RSAD2 (Radical S-Adenosyl Methionine Domain Containing 2) is a Protein Coding gene. Diseases associated with RSAD2 include Yellow Fever and Influenza. Among its related pathways are Interferon gamma signaling and Innate Immune System. GO annotations related to this gene include protein self-association and iron-sulfur cluster binding.

抗原:Synthesized peptide derived from the N-terminal region of human RSAD2.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: