

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey CBLN4 Polyclonal Antibody

多克隆抗体

产品货号: PA02689

第三版

描述: This gene encodes a member of a family of small secreted proteins containing C1Q domains. Members of this family are involved in regulation of neurexin signalling during synapse development. The mouse homolog of the protein encoded by this gene competes with netrin to bind to the deleted in colorectal cancer receptor. CBLN4 (Cerebellin 4 Precursor) is a Protein Coding gene. An important paralog of this gene is CBLN2.

抗原: Synthesized peptide derived from the C-terminal region of human Cerebellin 4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: