

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human FGR Polyclonal Antibody

多克隆抗体

产品货号: PA02696

第三版

**描述:** This gene is a member of the Src family of protein tyrosine kinases (PTKs). The encoded protein contains N-terminal sites for myristylation and palmitoylation, a PTK domain, and SH2 and SH3 domains which are involved in mediating protein-protein interactions with phosphotyrosine-containing and proline-rich motifs, respectively. The protein localizes to plasma membrane ruffles, and functions as a negative regulator of cell migration and adhesion triggered by the beta-2 integrin signal transduction pathway. Infection with Epstein-Barr virus results in the overexpression of this gene. Multiple alternatively spliced variants, encoding the same protein, have been identified.

**抗原:** Synthesized peptide derived from human c-Fgr around the non-phosphorylation site of Y412.

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Detected in neutrophils, monocytes and natural killer cells (at protein level). Detected in monocytes and large lymphocytes.