



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey CCZ1/CCZ1B Polyclonal Antibody

多克隆抗体

产品货号: PA02724

第三版

**描述:**CCZ1 (CCZ1 Homolog, Vacuolar Protein Trafficking And Biogenesis Associated) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and RAB GEFs exchange GTP for GDP on RABs. An important paralog of this gene is CCZ1B.

**抗原:**Synthesized peptide derived from the C-terminal region of human CGI-43

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**