

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CIB3 Polyclonal Antibody

多克隆抗体

产品货号: PA02809

第三版

描述: This gene product shares a high degree of sequence similarity with DNA-dependent protein kinase catalytic subunit-interacting protein 2 in human and mouse, and like them may bind the catalytic subunit of DNA-dependent protein kinases. The exact function of this gene is not known. Alternative splicing results in multiple transcript variants. CIB3 (Calcium And Integrin Binding Family Member 3) is a Protein Coding gene. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CIB2.

抗原: Synthesized peptide derived from the C-terminal region of human CIB3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: