

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Claudin 4 Polyclonal Antibody

多克隆抗体

产品货号: PA02835

第三版

描述: Claudin 4, also known as CLDN4, is a protein which in humans is encoded by the CLDN4 gene. It belongs to the group of claudins. This gene encodes an integral membrane protein, which belongs to the claudin family. The protein is a component of tight junction strands and may play a role in internal organ development and function during pre- and postnatal life. This gene is deleted in Williams-Beuren syndrome, a neurodevelopmental disorder affecting multiple systems.

抗原: Synthesized peptide derived from human Claudin-4 around the non-phosphorylation site of Y208.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: