
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CLCN4 Polyclonal Antibody

多克隆抗体

产品货号: PA02855

第三版

描述:The CLCN family of voltage-dependent chloride channel genes comprises nine members (CLCN1-7, Ka and Kb) which demonstrate quite diverse functional characteristics while sharing significant sequence homology. Chloride channel 4 has an evolutionary conserved CpG island and is conserved in both mouse and hamster. This gene is mapped in close proximity to APXL (Apical protein Xenopus laevis-like) and OA1 (Ocular albinism type I), which are both located on the human X chromosome at band p22.3. The physiological role of chloride channel 4 remains unknown but may contribute to the pathogenesis of neuronal disorders. Alternate splicing results in two transcript variants that encode different proteins. CLCN4 (Chloride Voltage-Gated Channel 4) is a Protein Coding gene. Diseases associated with CLCN4 include Mental Retardation, X-Linked 49 and Non-Syndromic X-Linked Intellectual Disability. Among its related pathways are Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds and Ion channel transport. GO annotations related to this gene include ion channel activity and antiporter activity. An important paralog of this gene is CLCN5.

抗原:Synthesized peptide derived from the Internal region of human CLC-4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Abundant in skeletal muscle and also detectable in brain and heart.