

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CLCNKA Polyclonal Antibody

多克隆抗体

产品货号: PA02867

第三版

**描述:** This gene is a member of the CLC family of voltage-gated chloride channels. The encoded protein is predicted to have 12 transmembrane domains, and requires a beta subunit called barttin to form a functional channel. It is thought to function in salt reabsorption in the kidney and potassium recycling in the inner ear. The gene is highly similar to CLCNKB, which is located 10 kb downstream from this gene. Multiple transcript variants encoding different isoforms have been found for this gene.

**抗原:** Synthesized peptide derived from the C-terminal region of human CLC-KA

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Expressed predominantly in the kidney. All nephron segments expressing BSND also express CLCNK proteins.