
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TP53BP1 Polyclonal Antibody

多克隆抗体

产品货号: PA02892

第三版

描述:TP53BP1 (tumor protein p53 binding protein 1) is a protein-coding gene. Go annotations related to this gene include methylated histone residue binding and p53 binding.Plays a key role in the response to DNA damage. May have a role in checkpoint signaling during mitosis.Enhances TP53-mediated transcriptional activation.

抗原:Synthesized peptide derived from human 53BP1 around the non-phosphorylation site of Ser6.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: