

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CMTM4 Polyclonal Antibody

多克隆抗体

产品货号: PA02918

第三版

描述: This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants encoding different isoforms have been identified. CMTM4 (CKLF Like MARVEL Transmembrane Domain Containing 4) is a Protein Coding gene. GO annotations related to this gene include cytokine activity. An important paralog of this gene is CMTM6.

抗原: Synthesized peptide derived from the Internal region of human CMTM4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in testis and prostate.