

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse VMA21 Polyclonal Antibody

多克隆抗体

产品货号: PA02921

第三版

描述: This gene encodes a chaperone for assembly of lysosomal vacuolar ATPase. Required for the assembly of the V0 complex of the vacuolar ATPase (V-ATPase) in the endoplasmic reticulum. Associates with the V0 complex of the vacuolar ATPase (V-ATPase). MEAX is a childhood-onset disease characterized by progressive vacuolation and atrophy of skeletal muscle. It is inherited in recessive fashion, affecting boys and sparing carrier females. Onset is in childhood, and patients exhibit weakness of the proximal muscles of the lower extremities, progressing slowly to involve other skeletal muscle groups over time.

抗原: Synthetic peptide of human VMA21

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: