

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC18A1 Polyclonal Antibody

多克隆抗体

产品货号: PA02924

第三版

描述:The vesicular monoamine transporter acts to accumulate cytosolic monoamines into vesicles, using the proton gradient maintained across the vesicular membrane. Its proper function is essential to the correct activity of the monoaminergic systems that have been implicated in several human neuropsychiatric disorders. The transporter is a site of action of important drugs, including reserpine and tetrabenazine.

抗原:Synthetic peptide of human SLC18A1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in chromaffin cells of the adrenal medulla (at protein level). Detected in peripheral sympathetic ganglia (at protein level). Found in some paracrine cells in stomach and duodenum (at protein level).1 Publication