



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZNF471 Polyclonal Antibody

多克隆抗体

产品货号: PA02934

第三版

**描述:** This protein belongs to the krueppel C2H2-type zinc-finger protein family. It may be involved in transcriptional regulation. The subcellular location of this protein is Nucleus. It contains 15 C2H2-type zinc fingers and one KRAB domain.

**抗原:** Synthetic peptide of human ZNF471

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**