

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ABCD1 Polyclonal Antibody

多克隆抗体

产品货号: PA02935

第三版

描述:The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the ALD subfamily, which is involved in peroxisomal import of fatty acids and/or fatty acyl-CoAs in the organelle. All known peroxisomal ABC transporters are half transporters which require a partner half transporter molecule to form a functional homodimeric or heterodimeric transporter. This peroxisomal membrane protein is likely involved in the peroxisomal transport or catabolism of very long chain fatty acids. Defects in this gene have been identified as the underlying cause of adrenoleukodystrophy, an X-chromosome recessively inherited demyelinating disorder of the nervous system. ABCD1 (ATP Binding Cassette Subfamily D Member 1) is a Protein Coding gene. Diseases associated with ABCD1 include Adrenoleukodystrophy and Addison's Disease. Among its related pathways are Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds and Peroxisome. GO annotations related to this gene include protein homodimerization activity and enzyme binding. An important paralog of this gene is ABCD2.

抗原:Synthesized peptide derived from the Internal region of human ABCD1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: