
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse UNC5CL Polyclonal Antibody

多克隆抗体

产品货号: PA02942

第三版

描述:This protein inhibits NF-kappa-B-dependent transcription by impairing NF-kappa-B binding to its targets. It interacts with p65/RELA and NFKB1. It has high levels expression in pancreas, liver and kidney. It contains one death domain and one ZU5 domain.

抗原:Synthetic peptide of human UNC5CL

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in pancreas, liver and kidney.