

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ABHD2 Polyclonal Antibody

多克隆抗体

产品货号: PA02955

第三版

描述:ABHD2 encodes a protein containing an alpha/beta hydrolase fold, which is a catalytic domain found in a wide range of enzymes. The encoded protein is an acylglycerol lipase that catalyzes the hydrolysis of endocannabinoid arachidonoylglycerol from the cell membrane. This leads to activation of the sperm calcium channel CatSper, which results in sperm activation. Alternative splicing of this gene results in two transcript variants encoding the same protein. ABHD2 (Abhydrolase Domain Containing 2) is a Protein Coding gene. GO annotations related to this gene include carboxylic ester hydrolase activity. An important paralog of this gene is ABHD3.

抗原:Synthesized peptide derived from the Internal region of human ABHD2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Present in sperm (at protein level).