

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse RCC1L Polyclonal Antibody

多克隆抗体

产品货号: PA02961

第三版

描述: This gene encodes a protein containing regulator of chromosome condensation 1-like repeats. The encoded protein may function as a guanine nucleotide exchange factor. This gene is located in a region of chromosome 7 that is deleted in Williams-Beuren syndrome, a multisystem developmental disorder. Alternative splicing results in multiple transcript variants. RCC1L (RCC1 Like) is a Protein Coding gene. Diseases associated with RCC1L include Williams-Beuren Syndrome. An important paralog of this gene is RPGR.

抗原: Synthetic peptide of human WBSCR16

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous.