

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CNT2 Polyclonal Antibody

多克隆抗体

产品货号: PA02965

第三版

描述:SLC28A2 (Solute Carrier Family 28 Member 2) is a Protein Coding gene. Among its related pathways are Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds and Thiopurine Pathway, Pharmacokinetics/Pharmacodynamics. GO annotations related to this gene include nucleoside binding and purine nucleoside transmembrane transporter activity. An important paralog of this gene is SLC28A1.

抗原:Synthesized peptide derived from the Internal region of human CNT2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in heart and skeletal muscle followed by liver, kidney, intestine, pancreas, placenta and brain. Weak expression in lung.