

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ACBD6 Polyclonal Antibody

多克隆抗体

产品货号: PA02975

第三版

描述:ACBD6 (Acyl-CoA Binding Domain Containing 6) is a Protein Coding gene. Among its related pathways are Fatty Acyl-CoA Biosynthesis and Metabolism. GO annotations related to this gene include lipid binding and fatty-acyl-CoA binding. An important paralog of this gene is OSTF1.

抗原:Synthesized peptide derived from the Internal region of human ACBD6

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in placenta and spleen (at protein level). Detected in placenta, umbilical cord blood, CD34-positive hematopoietic progenitor cells and bone marrow.