

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BVES Polyclonal Antibody

多克隆抗体

产品货号: PA02977

第三版

描述:This gene encodes a member of the POP family of proteins containing three putative transmembrane domains. This gene is expressed in cardiac and skeletal muscle and may play an important role in development of these tissues. The mouse ortholog may be involved in the regeneration of adult skeletal muscle and may act as a cell adhesion molecule in coronary vasculogenesis. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Dec 2010]

抗原:Recombinant protein of human BVES

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in epithelial cells (at protein level). Expressed in fetal and adult heart and skeletal muscle.