

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human COL4A6 Polyclonal Antibody

多克隆抗体

产品货号: PA02993

第三版

**描述:** This gene encodes one of the six subunits of type IV collagen, the major structural component of basement membranes. Like the other members of the type IV collagen gene family, this gene is organized in a head-to-head conformation with another type IV collagen gene, alpha 5 type IV collagen, so that the gene pair shares a common promoter. Deletions in the alpha 5 gene that extend into the alpha 6 gene result in diffuse leiomyomatosis accompanying the X-linked Alport syndrome caused by the deletion in the alpha 5 gene. Alternative splicing results in multiple transcript variants encoding different isoforms. COL4A6 (Collagen Type IV Alpha 6 Chain) is a Protein Coding gene. Diseases associated with COL4A6 include Deafness, X-Linked 6 and X-Linked Non-Syndromic Sensorineural Deafness Type Dfn. Among its related pathways are Collagen chain trimerization and Phospholipase-C Pathway. GO annotations related to this gene include structural molecule activity and extracellular matrix structural constituent. An important paralog of this gene is COL4A2.

**抗原:** Synthesized peptide derived from the C-terminal region of human COL4A6

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**