

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ACOT4 Polyclonal Antibody

多克隆抗体

产品货号: PA02997

第三版

描述:ACOT4 (Acyl-CoA Thioesterase 4) is a Protein Coding gene. Among its related pathways are Fatty acid elongation and Fatty Acyl-CoA Biosynthesis. GO annotations related to this gene include receptor binding and palmitoyl-CoA hydrolase activity. An important paralog of this gene is ACOT1.

抗原:Synthesized peptide derived from the C-terminal region of human ACOT4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Strongest expression in liver and kidney and weaker expression in placenta, heart, and muscle.