

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Connexin 25 Polyclonal Antibody

多克隆抗体

产品货号: PA03004

第三版

描述:Connexins, such as GJB7, are involved in the formation of gap junctions, intercellular conduits that directly connect the cytoplasm of contacting cells. Each gap junction channel is formed by docking of 2 hemichannels, each of which contains 6 connexin subunits. GJB7 (Gap Junction Protein Beta 7) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and Gap junction trafficking. An important paralog of this gene is GJB5.

抗原:Synthesized peptide derived from the Internal region of human Connexin 25

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Weakly expressed in placenta.