

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ADIPOQ Polyclonal Antibody

多克隆抗体

产品货号: PA03005

第三版

描述:Rat gAcrp30 (amp-1) is a 16.6 kDa protein that corresponds to the C-terminal globular domain of adiponectin (amp-1). gAcrp30/Adipolean is expressed and secreted exclusively by adipose tissue. This protein possesses unique pharmacological properties for regulation of body weight and lipid metabolism. It is a key component in controlling energy homeostasis. Recently, gAcrp30/Adipolean has been shown to inhibit the proliferation of myelomonocytic progenitors. The Recombinant rat gAcrp30/Adipolean is a 16.6 kDa protein consisting of 145 amino acid residues

抗原:Synthesized peptide derived from the Internal region of human Acrp30.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Synthesized exclusively by adipocytes and secreted into plasma.