



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse RPA1 Polyclonal Antibody

多克隆抗体

产品货号: PA03008

第三版

**描述:**RPA1 (replication protein A1, 70kDa) is a protein-coding gene. Diseases associated with RPA1 include xeroderma pigmentosum, group a, and werner syndrome. GO annotations related to this gene include chromatin binding and single-stranded DNA binding.

**抗原:**Recombinant protein of human RPA1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**