

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLC27A5 Polyclonal Antibody

多克隆抗体

产品货号: PA03015

第三版

描述:The protein encoded by this gene is an isozyme of very long-chain acyl-CoA synthetase (VLCS). It is capable of activating very long-chain fatty-acids containing 24- and 26-carbons. It is expressed in liver and associated with endoplasmic reticulum but not with peroxisomes. Its primary role is in fatty acid elongation or complex lipid synthesis rather than in degradation. This gene has a mouse ortholog.

抗原:Synthesized peptide derived from the Internal region of human ACSVL6

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Predominantly expressed in liver.