

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Connexin 47 Polyclonal Antibody

多克隆抗体

产品货号: PA03018

第三版

描述:One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell. May play a role in myelination in central and peripheral nervous systems.

抗原:Synthesized peptide derived from the N-terminal region of human Connexin 47

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in central nervous system, in sciatic nerve and sural nerve. Also detected in skeletal muscles.