

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CPA5 Polyclonal Antibody

多克隆抗体

产品货号: PA03060

第三版

描述: Carboxypeptidases have functions ranging from digestion of food to selective biosynthesis of neuroendocrine peptides. Members of the A/B subfamily of carboxypeptidases, such as CPA5, contain an approximately 90-amino acid pro region that assists in the folding of the active carboxypeptidase domain. Cleavage of the pro region activates the enzyme. CPA5 (Carboxypeptidase A5) is a Protein Coding gene. Diseases associated with CPA5 include Amebiasis and Silver-Russell Syndrome. GO annotations related to this gene include peptidase activity and carboxypeptidase activity. An important paralog of this gene is CPA1.

抗原: Synthesized peptide derived from the Internal region of human CPA5

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expression is very low or not detectable.