

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CPM Polyclonal Antibody

多克隆抗体

产品货号: PA03068

第三版

描述:Specifically removes C-terminal basic residues (Arg or Lys) from peptides and proteins. It is believed to play important roles in the control of peptide hormone and growth factor activity at the cell surface, and in the membrane-localized degradation of extracellular proteins.

抗原:Synthesized peptide derived from the Internal region of human CPM

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: