

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AIG1 Polyclonal Antibody

多克隆抗体

产品货号: PA03074

第三版

描述:AIG1 (Androgen Induced 1) is a Protein Coding gene. An important paralog of this gene is ADTRP.

抗原:Synthesized peptide derived from the C-terminal region of human AIG-1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in heart, ovary, testis, liver, and kidney, at lower levels in spleen, prostate, brain, skeletal muscle, pancreas, small intestine and colon, and undetected in peripheral blood leukocytes, thymus, lung and placenta. AIG1 expression is higher in hair follicles from males than from females.