

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat AKAP10 Polyclonal Antibody

多克隆抗体

产品货号: PA03085

第三版

描述: This gene encodes a member of the A-kinase anchor protein family. A-kinase anchor proteins bind to the regulatory subunits of protein kinase A (PKA) and confine the holoenzyme to discrete locations within the cell. The encoded protein is localized to mitochondria and interacts with both the type I and type II regulatory subunits of PKA. Polymorphisms in this gene may be associated with increased risk of arrhythmias and sudden cardiac death. AKAP10 (A-Kinase Anchoring Protein 10) is a Protein Coding gene. Diseases associated with AKAP10 include Cardiac Conduction Defect and Sudden Cardiac Death. Among its related pathways are Activation of cAMP-Dependent PKA and Factors involved in megakaryocyte development and platelet production.

抗原: Synthesized peptide derived from the N-terminal region of human AKAP 10

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: