

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CSP Polyclonal Antibody

多克隆抗体

产品货号: PA03143

第三版

描述: This gene is a member of the J protein family. J proteins function in many cellular processes by regulating the ATPase activity of 70 kDa heat shock proteins. The encoded protein plays a role in membrane trafficking and protein folding, and has been shown to have anti-neurodegenerative properties. The encoded protein is known to play a role in cystic fibrosis and Huntington's disease. A pseudogene of this gene is located on the short arm of chromosome 8.

抗原: Synthesized peptide derived from the C-terminal region of human CSP

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed in pancreas, kidney, skeletal muscle, liver, lung, placenta, brain and heart.