

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat/Monkey Cyclin A Polyclonal Antibody

多克隆抗体

产品货号: PA03198

---

第三版

**描述:**Essential for the control of the cell cycle at the G1/S (start) and the G2/M (mitosis) transitions.

**抗原:**Synthesized peptide derived from the Internal region of human Cyclin A

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Very high levels in testis and very low levels in brain. Also found in myeloid leukemia cell lines.