
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat HNRNP K Polyclonal Antibody

多克隆抗体

产品货号: PA03206

第三版

描述:One of the major pre-mRNA-binding proteins. Binds tenaciously to poly(C) sequences. Likely to play a role in the nuclear metabolism of hnRNAs, particularly for pre-mRNAs that contain cytidine-rich sequences. Can also bind poly(C) single-stranded DNA.

抗原:Recombinant protein of human HNRNPK

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: