

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ALDH3B1 Polyclonal Antibody

多克隆抗体

产品货号: PA03239

第三版

**描述:**The aldehyde dehydrogenases are a family of isozymes that may play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. This particular gene spans about 20 kb of genomic DNA and is composed of 9 coding exons. The gene is highly expressed in kidney and lung. The functional significance of this gene as well as the cellular localization of its product are presently unknown. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

**抗原:**Synthesized peptide derived from the Internal region of human ALDH3B1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$  g/ml.

**特异性:**Highest expression in kidney and lung.