

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ALX3 Polyclonal Antibody

多克隆抗体

产品货号: PA03258

第三版

描述:This gene encodes a nuclear protein with a homeobox DNA-binding domain that functions as a transcriptional regulator involved in cell-type differentiation and development. Preferential methylation of this gene's promoter is associated with advanced-stage neuroblastoma tumors.ALX3 (ALX Homeobox 3) is a Protein Coding gene. Diseases associated with ALX3 include Frontonasal Dysplasia 1 and Basal Encephalocele. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is ALX4.

抗原:Synthesized peptide derived from the Internal region of human ALX3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: