

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat CD83 Polyclonal Antibody

多克隆抗体

产品货号: PA03272

第三版

**描述:**The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene. CD83 (CD83 Molecule) is a Protein Coding gene. Diseases associated with CD83 include Exanthema Subitum and Specific Bursitis Often Of Occupational Origin. Among its related pathways are T cell receptor signaling pathway and B Cell Development Pathways.

**抗原:**Recombinant protein of human CD83

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed by activated lymphocytes, Langerhans cells and interdigitating reticulum cells.