

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANKRD30A Polyclonal Antibody

多克隆抗体

产品货号: PA03279

第三版

描述: This gene encodes a DNA-binding transcription factor that is uniquely expressed in mammary epithelium and the testis. Altered expression levels have been associated with breast cancer progression. ANKRD30A (Ankyrin Repeat Domain 30A) is a Protein Coding gene. Diseases associated with ANKRD30A include Breast Cancer. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding. An important paralog of this gene is ANKRD30B.

抗原: Synthesized peptide derived from the C-terminal region of human ANKRD30A

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Mainly expressed in breast and testis. A very faint signal is detected in placenta. Also expressed in many breast cancer cells.