

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NUSAP1 Polyclonal Antibody

多克隆抗体

产品货号: PA03280

第三版

描述:NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization. May associate with chromosomes and promote the organization of mitotic spindle microtubules around them.

抗原:Synthesized peptide derived from the C-terminal region of human ANKT

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: