
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat NCF4 Polyclonal Antibody

多克隆抗体

产品货号: PA03291

第三版

描述:The protein encoded by this gene is a cytosolic regulatory component of the superoxide-producing phagocyte NADPH-oxidase, a multicomponent enzyme system important for host defense. This protein is preferentially expressed in cells of myeloid lineage. It interacts primarily with neutrophil cytosolic factor 2 (NCF2/p67-phox) to form a complex with neutrophil cytosolic factor 1 (NCF1/p47-phox), which further interacts with the small G protein RAC1 and translocates to the membrane upon cell stimulation. This complex then activates flavocytochrome b, the membrane-integrated catalytic core of the enzyme system. The PX domain of this protein can bind phospholipid products of the PI(3) kinase, which suggests its role in PI(3) kinase-mediated signaling events. The phosphorylation of this protein was found to negatively regulate the enzyme activity. Alternatively spliced transcript variants encoding distinct isoforms have been observed.NCF4 (Neutrophil Cytosolic Factor 4) is a Protein Coding gene. Diseases associated with NCF4 include Granulomatous Disease, Chronic, Autosomal Recessive, Cytochrome B-Positive, Type Iii and Chronic Granulomatous Disease. Among its related pathways are Toxoplasmosis and Blood-Brain Barrier and Immune Cell Transmigration: VCAM-1/CD106 Signaling Pathways. GO annotations related to this gene include protein dimerization activity and phosphatidylinositol-3-phosphate binding.

抗原:Recombinant protein of human NCF4

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Expression is restricted to hematopoietic cells.