



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat LETMD1 Polyclonal Antibody

多克隆抗体

产品货号: PA03306

第三版

描述:This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene.

抗原:Recombinant protein of human LETMD1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Kidney, liver, skeletal muscle, heart and brain. Overexpressed in various tumors including leukemia, lymphoma, and carcinomas of the breast, kidney, ovary, stomach, colon and uterine cervix.1 Publication