

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CYP26A1 Polyclonal Antibody

多克隆抗体

产品货号: PA03373

第三版

描述: This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This endoplasmic reticulum protein acts on retinoids, including all-trans-retinoic acid (RA), with both 4-hydroxylation and 18-hydroxylation activities. This enzyme regulates the cellular level of retinoic acid which is involved in regulation of gene expression in both embryonic and adult tissues. Two alternatively spliced transcript variants of this gene, which encode the distinct isoforms, have been reported. CYP26A1 (Cytochrome P450 Family 26 Subfamily A Member 1) is a Protein Coding gene. Diseases associated with CYP26A1 include Caudal Regression Syndrome and Mixed Germ Cell Cancer. Among its related pathways are Signaling by Retinoic Acid and Cytochrome P450 - arranged by substrate type. GO annotations related to this gene include iron ion binding and oxidoreductase activity, acting on paired donors, with incorporation or reduction of molecular oxygen. An important paralog of this gene is CYP26C1.

抗原: Synthesized peptide derived from the Internal region of human CYP26A1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highest levels in adult liver, heart, pituitary gland, adrenal gland, placenta and regions of the brain.