

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse /Rat EIF2 alpha Polyclonal Antibody

多克隆抗体

产品货号: PA03382

第三版

描述: Functions in the early steps of protein synthesis of a small number of specific mRNAs. Acts by directing the binding of methionyl-tRNA_i to 40S ribosomal subunits. In contrast to the eIF-2 complex, it binds methionyl-tRNA_i to 40 S subunits in a codon-dependent manner, whereas the eIF-2 complex binds methionyl-tRNA_i to 40 S subunits in a GTP-dependent manner. May act by impinging the expression of specific proteins.

抗原: Synthetic peptide of human EIF2S1 (Ab-51)

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: