

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human /Mouse FAM129B Polyclonal Antibody

多克隆抗体

产品货号: PA03396

第三版

描述:FAM129B (Family With Sequence Similarity 129 Member B) is a Protein Coding gene. Among its related pathways are MAP Kinase Signaling. An important paralog of this gene is FAM129A.

抗原:Synthetic peptide of human FAM129B (Ab-712)

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: