

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Monkey CYP2J2 Polyclonal Antibody

多克隆抗体

产品货号: PA03399

第三版

描述: This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and is thought to be the predominant enzyme responsible for epoxidation of endogenous arachidonic acid in cardiac tissue. Multiple transcript variants have been found for this gene. CYP2J2 (Cytochrome P450 Family 2 Subfamily J Member 2) is a Protein Coding gene. Diseases associated with CYP2J2 include Clopidogrel Resistance. Among its related pathways are Cytochrome P450 - arranged by substrate type and Arachidonate Epoxygenase / Epoxide Hydrolase. GO annotations related to this gene include iron ion binding and oxidoreductase activity, acting on paired donors, with incorporation or reduction of molecular oxygen. An important paralog of this gene is CYP2R1.

抗原: Synthesized peptide derived from the Internal region of human CYP2J2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Highly expressed in heart, present at lower levels in liver, ileum, jejunum, colon, and kidney.